

Thyroid Ultrasound - Image Protocol

1. Transverse midline isthmus
2. Transverse midline isthmus with AP measurement
3. Transverse RIGHT lobe upper pole
4. Transverse RIGHT lobe mid
5. Transverse RIGHT lobe mid with width measurement
6. Transverse RIGHT lobe mid with color Doppler
7. Transverse RIGHT lobe lower pole
8. Long RIGHT lobe lateral
9. Long RIGHT lobe mid
10. Long RIGHT lobe with length/AP measurements
11. Long RIGHT lobe medial
12. Transverse LEFT lobe upper pole
13. Transverse LEFT lobe mid
14. Transverse LEFT lobe mid with measurement
15. Transverse LEFT lobe mid with color Doppler
16. Transverse LEFT lobe lower pole
17. Long LEFT lobe lateral
18. Long LEFT lobe mid
19. Long LEFT lobe mid with measurement
20. Long LEFT lobe medial